

Screen Name: 1.0 Home Page / Login

Screen Design Layout: Note: The screen design layout does not represent the Graphic Design nor the font selection and size.

TOP MENU	Home	CareHelper Guide	About Us	Contact Us	Site Map
----------	------	------------------	----------	------------	----------

< Error Messages Screen Area: Text to be displayed in red>

Welcome to CareHelper.com

If you are already a CareHelper User, welcome back!

Please enter your CareGiver User Name and Password.

* CareGiver User Name

* Password

☒ (Save my password)

* Indicates required field entry

Legend: * Indicates required field entry; defines field entry; represents screen buttons;

Fig. 1

Screen Name: 2.0 New CareGiver Registration

Screen Design Layout: Note: The screen design layout does not represent the Graphic Design nor the font selection and size.

TOP MENU	Home	Benefits of Membership	Join Now	TBD	Contact Us
----------	------	------------------------	----------	-----	------------

< Error Messages Screen Area: Text to be displayed in red>

New CareGiver Registration

Just a few questions, and you are registered!

* CareGiver First Name

* CareGiver Last Name

* CareGiver Birth Date

* CareGiver Email Address

Please check your Email Address! The confirmation of your registration will be sent to this address.

Choose your CareGiver User Name and Password. This will identify you in CareHelper.com

* CareGiver User Name

☐ (Save my User Name, so I don't have to remember it each time I login)

Enter your password and re-type it to confirm. * Indicates required field entry

* Password Your password must be at least 4 characters in length

* Confirm Password

☐ (Save my password, so I don't have to remember it each time I login)

Just one more thing,

As a registered CareGiver of CareHelper.com, you agree to abide by our Terms of Use. Please read the following Terms of Use before submitting your information.

* Indicates required field entry

Legend: * Indicates required field entry; defines field entry; represents screen buttons;

Screen Design Layout:

Note: The screen design layout does not represent the Graphic Design nor the font selection and size.

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us
----------	------	--------------------	------------------	------------------	------------

< Error Messages Screen Area: Text to be displayed in red >

CareGiver Page

Welcome, <caregiver.firstname>, to your personalized Home Page.

If you are not Janice Smith, press here!

Alerts

<patient.nickname>

<event.start> <event.description>

<event.start> <event.description>

<patient.nickname>

<event.start> <event.description>

<event.start> <event.description>

Select or Add a Care Recipient,

Mom	May Smith	Active
Dad	Joseph P. Smith, Sr.	Inactive

<patient.nickname> <patient.fullname> <patient.status>

Add New Care Recipient

[Select] [Change Status] [Delete]

Shopping Manager

<caregiver.firstname>'s List

<shoplist.item_desc>

<shoplist.item_desc>

<patient.nickname>'s List

<shoplist.item_desc>

<shoplist.item_desc>

<patient.nickname>'s List

<shoplist.item_desc>

<shoplist.item_desc>

Legend: * Indicates required field entry;
 <table.field name> defines field entry;
 buttons represents screen buttons;

Screen Name: 3.1 Caregiver Personal Information

Screen Layout Design:

TOP MENU

Home

Janice's Organizer

Shopping Manager

CareHelper Guide

Contact Us

< error messages screen area >

<caregiver.nickname>'s Personal Information



Please complete the following information and bring the information current:

*First Name	<input type="text" value="<caregiver.firstname >"/>	*Last Name	<input type="text" value="<caregiver.lastname>"/>
Nickname	<input type="text" value="<caregiver.nickname >"/>	*Birth Date	<input type="text" value="<caregiver.birthdate> (mm/dd/yyyy)"/>
Gender	<input type="text" value="<caregiver.gender >"/>		
Address	<input type="text" value="<caregiver.address1 >"/> <input type="text" value="<caregiver.address2 >"/>		
City	<input type="text" value="<caregiver.city >"/>		
State	<input type="text" value="<caregiver.state >"/>	Zip Code	<input type="text" value="<caregiver.zip >"/>
*Day Phone	<input type="text" value="<caregiver.homephone>"/>	Mobile Phone	<input type="text" value="<caregiver.mobile >"/>
*Evening Phone	<input type="text" value="<caregiver.workphone >"/>	Pager No.	<input type="text" value="<caregiver.pagerno >"/>
Fax Phone	<input type="text" value="<caregiver.fax >"/>	Email	<input type="text" value="<caregiver.email >"/>
Primary Language	<input type="text" value="<caregiver.language1 >"/>	Secondary Language	<input type="text" value="<caregiver.language2>"/>

Notes:

* Indicates required field entry

[Save and Continue]


[Finish Later]

[Do not Save]

< error messages screen area >

Add New Care Recipient

PROCESS STEPS



Welcome, <caregiver.firstname>,

Please complete the following information about your Care Recipient:

* First Name

<patient.firstname>

* Last Name

<patient.lastname>


* Nickname

<patient.nickname> (eg. Mom, Dad, Nana, Patti, etc.)

* Birth Date

<patient.birthdate > (mm/dd/yyyy)

* Gender

<patient.gender > 

[Add New Care Recipient](#)

[Do not add, Return to Care Giver Page](#)

* Indicates required field entry

Fig. 5

Screen Name: 4.1 Care Recipient Personal Information

Screen Layout Design:

TOP MENU

Home

Janice's Organizer

Shopping Manager

CareHelper Guide

Contact Us

< error messages screen area >

<patient.nickname>'s Personal Information

PROCESS STEPS

Add New Care Recipient
Personal Information
Emergency Contact Information
Medical Information
Monitor Devices
Medications
Supplies
Pharmacy / Store
Reports

Please complete the following information about <patient.nickname> and bring the information current:

*First Name	<patient.firstname >	*Last Name	<patient.lastname>
*Nickname	<patient.nickname >	*Birth Date	<patient.birthdate> (mm/dd/yyyy)
*Gender	<patient.gender >↓	SSN	<patient.ssn > (123-45-6789)
Address	<patient.address1 >		
	<patient.address2 >		
City	<patient.city >		
State	<patient.state >↓	Zip Code	<patient.zip >
Day Phone	<patient.phone1 >	Evening Phone	<patient.phone2 >
Mobile Phone	<patient.mobile >		
Pager No.	<patient.pagermo >	Notification Email	<patient.email >
Primary Language	<patient.language1 >↓	Secondary Language	<patient.language2>↓

Notes:

* Indicates required field entry

[Save and Continue]

[Finish Later]

[Do not Save]

Screen Layout Design:

TOP MENU

Home

Janice's Organizer

Shopping Manager


CareHelper Guide

Contact Us

< error messages screen area >

<patient.nickname>'s Emergency Contact Information

PROCESS STEPS



Please select an Emergency Contact and complete the following information:

Select or Add an Emergency Contact,

Primary Contact Order	Name	Relationship To Care Recipient
1)	Tom Comings	Son
2)	Judith Comings	Neighbor
New Contact		

* Contact Order

<contacts.order> , in case of an emergency

First Name

<contacts.firstname>

Last Name

<contacts.lastname>

Nickname

<contacts.nickname>

Birth Date

<contacts.birthdate> (mm/dd/yyyy)

* Relationship

<contacts.relationship>

(i.e. Relationship to Care Recipient)

* Evening Phone

<contacts.home>

* Day Phone

<contacts.work>

Mobile Phone

<contacts.mobile>

Pager No.

<contacts.pagerno>

Email

<contacts.email>

Primary

<contacts.language1>

Secondary

<contact.language2>

Language

<contacts.language1>

Notes: (eg. Primary Contact has security code for home)

* Indicates required field entry

Save and Continue

Finish Later

Do not Save

Screen Name: 4.3.3 Care Recipient Physician Information

Screen Layout Design:

TOP MENU Home Janice's Organizer Shopping Manager CareHelper Guide Contact Us

< error messages screen area >

<patient.nickname>'s Medical - Physician Information

Please select a Physician and complete the following information:

PROCESS STEPS

Add New Care Recipient

Emergency Contact Information

Medical Information

ADL

Medical History

Physician

Monitor Devices

Medications

Supplies

Pharmacy / Store

Reports

Select or Add a Physician,

Name Specialty Status

Tom Comings, M.D.	Geriatrics	Active
Sally Jones, M.D.	Geriatrics	Inactive
Add New Physician		

Status <physician.status>

First Name <physician.firstname>

Last Name <physician.lastname>

Office Contact <physician.contact>

Office Phone <physician.office>

Fax Number <physician.fax>

Address <physician.address1>

City <physician.address2>

State <physician.city>

Zip Code <physician.zip>

Specialty 1 <physician.specialty1>

Other: <physician.specialty1>

Notes:

[Save and Continue]


[Finish Later]

[Do Not Save]

Screen Design Layout: Shopping Mgr

< error messages screen area >

Shopping Manager

Suggestions and Reminders				
[SORT]	[ADD ALL]	<u>Item</u>	<u>Person</u>	Delete 
		11/23/1999 <item_description	> <patient.nickname>	<input checked="" type="checkbox"/> Add Now
		11/19/1999 <item_description	> <patient.nickname>	<input checked="" type="checkbox"/> Add Now
		11/18/1999 <item_description	> <patient.nickname>	<input type="checkbox"/> Add Now
		11/18/1999 <item_description	> <patient.nickname>	<input type="checkbox"/> Add Now

Suggestions and Reminders

Prior Shopping Items

[SORT]	[ADD ALL]		<u>Del.</u>
<cargiver.firstname>'s List			
Medications			
<item_description	<input type="checkbox"/>	<input type="button" value="Add Now"/>	<input type="checkbox"/>
<item_description	<input type="checkbox"/>	<input type="button" value="Add Now"/>	<input type="checkbox"/>
Supplies			
<item_description	<input checked="" type="checkbox"/>	<input type="button" value="Add Now"/>	<input checked="" type="checkbox"/>
Other			
<item_description	<input type="checkbox"/>	<input type="button" value="Add Now"/>	<input type="checkbox"/>
<item_description	<input type="checkbox"/>	<input type="button" value="Add Now"/>	<input type="checkbox"/>
<patient.nickname>'s List			
Medications			
<item_description	<input type="checkbox"/>	<input type="button" value="Add Now"/>	<input type="checkbox"/>
<item_description	<input type="checkbox"/>	<input type="button" value="Add Now"/>	<input type="checkbox"/>
Supplies			
<item_description	<input checked="" type="checkbox"/>	<input type="button" value="Add Now"/>	<input checked="" type="checkbox"/>

Prior Shopping Items

<u><patient.nickname></u>		
<u>Qty</u>	<u>Item</u>	<u>Delete</u>
001 ↓	<shoplist.item desc	<input type="checkbox"/>
001 ↓	<shoplist.item desc	<input type="checkbox"/>
001 ↓	<shoplist.item desc	<input type="checkbox"/>
001 ↓	<shoplist.item desc	<input type="checkbox"/>
001 ↓	<shoplist.item desc	<input type="checkbox"/>
001 ↓	<shoplist.item desc	<input type="checkbox"/>
001 ↓	<shoplist.item desc	<input type="checkbox"/>

<u><patient.nickname></u>		
<u>Qty</u>	<u>Item</u>	<u>Delete</u>
001 ↓	<shoplist.item desc	<input type="checkbox"/>
001 ↓	<shoplist.item desc	<input type="checkbox"/>

<patient.nickname>'s List

Medications				
<item_description	>	<cg.location> ➡	Add Now	<input type="checkbox"/>
<item_description	>	<cg.location> ➡	Add Now	<input type="checkbox"/>
Supplies				
<item_description	>	<cg.location> ➡	Add Now	<input checked="" type="checkbox"/>

SAVE FOR LATER

EXIT

என்

Screen Name: 4.7 Monitor Device Set Up

Screen design Layout:

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us
----------	------	--------------------	------------------	------------------	------------

< error messages screen area >

<patient.nickname>'s Monitor Device Set Up

Select or Add a Device,

Device	Collection Method
Pill Bottle	Auto
Blood Pressure	na
Blood Glucose	Auto
Body Weigh	Manual
Dose Pager	na
New Device	

Fig. 10

Screen Name: 4.7.5 Pill Bottle Device Set Up

Screen design Layout:

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us
<p align="center">< error messages screen area ></p> <p align="center"><patient.nickname>'s Pill Bottle Device Set Up</p>					
Pill Bottle Set Up Instructions Step 1: Step 2: Step 3:	<p align="center">Select a Medication for Pill Bottle Set Up</p>				
	Medications	Dispense Method	Pill Bottle No.	Status	
	<patient.meds>	Pill Bottle	101	Active	
	<patient.meds>	Manual	na	Inactive	
	<patient.meds>	Pill Bottle	<pbottle.assignno>	<pbottle.status>	
<p align="center">New Device</p>					
<p align="center"><patient.meds></p>					
Dispense Method		Pill Bottle	Assigned Number	<pbottle.assignno>	
Status		<pbottle.status>			
Pill Bottle Begin Date		<pbottle.begindate>	Pill Bottle Stop Date		
Bottle Fill Days Supply		<dayssupply>	(i.e. number of pills placed in bottle for each fill)		
Dosage		<medsdosage>	Units	<medsunits>	
<p align="center">Times per day <medsinterval> or O As Needed</p>					
Set Daily Dosage Time:		12:15 AM	Alarm? • Yes ○ No		
Set Daily Dosage Time:		12:15 AM	Alarm? • Yes ○ No		
[Setup Pill Bottle]		[Finish Later]	[Refill]	[Do Not Save]	

Italic fields are display only (from medications set up screen)

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us										
<p align="center">< error messages screen area ></p> <h2 align="center"><patient.nickname>'s Monitor Device Set Up</h2>															
TBD	<p align="center">Select or Add Device,</p>														
...	<table border="1"> <tr> <td>Device</td> <td>Collection Method</td> </tr> <tr> <td>Weight/Temp Measure</td> <td>Manual</td> </tr> <tr> <td>Blood Pressure</td> <td>na</td> </tr> <tr> <td>Blood Glucose</td> <td>Auto</td> </tr> <tr> <td>New Device</td> <td></td> </tr> </table>					Device	Collection Method	Weight/Temp Measure	Manual	Blood Pressure	na	Blood Glucose	Auto	New Device	
Device	Collection Method														
Weight/Temp Measure	Manual														
Blood Pressure	na														
Blood Glucose	Auto														
New Device															
<h3 align="center"><u>Weight/Temp</u></h3> <p>Please select the date for <patient.nickname>'s Weight/Temp Set Up and complete the information:</p> <p>Weight/Temp <Date> <weight> [Add New Weight/Temp Measure]</p> <p>Please complete the following information:</p> <table border="1"> <tr> <td>Measure Method</td> <td>Manual</td> </tr> <tr> <td>Date</td> <td><Date></td> </tr> <tr> <td>Weight Measure</td> <td><wmeasure> Range: High <weighthigh> Low <weightlow></td> </tr> <tr> <td>Temp Measure</td> <td><tmeasure> Range: High <temphigh> Low <templow></td> </tr> </table> <p align="right"> [Save and Continue] [Finish Later] [Do Not Save] </p>						Measure Method	Manual	Date	<Date>	Weight Measure	<wmeasure> Range: High <weighthigh> Low <weightlow>	Temp Measure	<tmeasure> Range: High <temphigh> Low <templow>		
Measure Method	Manual														
Date	<Date>														
Weight Measure	<wmeasure> Range: High <weighthigh> Low <weightlow>														
Temp Measure	<tmeasure> Range: High <temphigh> Low <templow>														

Fig. 12

CareHelper.com Screen Flow




Fig. 13

Legend: CG - CareGiver; CR - Care Recipient;

Screen Name: NS - 1.1 Test Camera (resides on Nurse Station)

Step by Step Instructions

1. Adjust camera hardware
2. Preview video Adjust Volume, Contrast and Brightness
3. Preview Video
4. Save Setting

Test Camera

Fig. 14

01:30:00 Count Down 5 4 3 2 1




[PREVIEW]

[RECORD]


Microphone Level



[Volume ↑ ↓]



[Contrast ↑ ↓] [Brightness ↑ ↓]



[SAVE SETTING]

[EXIT]

[RESET]

Overview:

1. The purpose of this screen is to allow the Nurse to be able to view themselves and adjust the hardware and controls (Volume, Contrast and Brightness, only)
2. Test Camera Screen resides on the Nurse Station along with the Camera software.
3. The video is recorded as Mirror Image.
4. The camera is always on, see viewing the image starts immediately when first entering the screen.
5. [RECORD] does a count down from 5, before it starts recording
6. [PREVIEW] previews the last test recorded video.
7. Control settings show the current and recommended range of volume, contrast and brightness.
8. Above the Volume controls is a microphone recording feedback bar, which indicates the microphone recording level.
9. Note: The technical video settings that impact processing, quality, performance, compression and sending the video is to be accessed by a separate System Administrator screen only.

Screen Name: NS-1.2 Select CareGiver and Care Recipient (resides on web server)


Fig. 15

Select CareGiver		Select Care Recipients	
<caregiver.last>, <caregiver.first>	<Patient nickname> <patient.first> <patient.last>		
Carroll, Pat	Mom	John	Smith
LeMore, Nancy			
McCammon, Beth			
Otters, Joan			
Smith, Stephanie			
Thompson, Jane			
Victor, Andre			
Vace, Kenneth			
Young, Lee			
[Continue]		[EXIT]	

Overview:

1. Purpose of screen is to select who we are making the presentation for: CareGiver / Care Recipient
2. All of the CareGivers previously assigned to a nurse will be displayed in the CareGiver List. Default to highlighting the first CareGiver in the list.
3. When a CareGiver is highlighted, all of the Care Recipients associated with the CareGiver will be displayed in the Care Recipient window.
 - default highlighting the first Care Recipient listed.
4. This screen is based on the web-server, and accessed by the nurse through the Nurse Base Station. This interface is has seamless as possible to the nurse.
5. After selecting the CG, CR and the[Continue] button is pressed, the selected information required to complete the Presentation is returned to the Nurse Station, to complete the video.
6. EXIT exits the screen and returns to the NS Main Screen.

Screen Name: NS - 1.3 Select Reports for Presentation (resides on web server)

CareGiver: Eric Brown Care Recipient: Mary Smith		Presentation Date: mm/dd/yyyy Nurse: Ms. Kathy Kindess																
Select Reports to include in Presentation: [View Prior Presentation] <table><thead><tr><th>Include</th><th>Reports</th><th>Description</th></tr></thead><tbody><tr><td><input checked="" type="checkbox"/></td><td>Report Template 1</td><td>Report Template 1 Descr</td></tr><tr><td><input type="checkbox"/></td><td>Report Template 2</td><td></td></tr><tr><td><input checked="" type="checkbox"/></td><td>Report Template 3</td><td></td></tr><tr><td><input checked="" type="checkbox"/></td><td>Alerts</td><td></td></tr></tbody></table>		Include	Reports	Description	<input checked="" type="checkbox"/>	Report Template 1	Report Template 1 Descr	<input type="checkbox"/>	Report Template 2		<input checked="" type="checkbox"/>	Report Template 3		<input checked="" type="checkbox"/>	Alerts		Presentation Summary Bullet Points: <ul style="list-style-type: none">Blood pressure in normal rangeetc.	
Include	Reports	Description																
<input checked="" type="checkbox"/>	Report Template 1	Report Template 1 Descr																
<input type="checkbox"/>	Report Template 2																	
<input checked="" type="checkbox"/>	Report Template 3																	
<input checked="" type="checkbox"/>	Alerts																	
View Selected Report: 		Enter Clinical Notes for Selected Report: Patient periodically is outside the recommended range, and needs to be monitored closely during the holidays.																
<div>[<< BACK] [NEXT >>] [EXIT] [RESET]</div>																		

Overview:

- Purpose of screen is to select the reports to be included in the presentation. In addition, clinical notes, that remain confidential to the clinician may be entered for each report
- "1.3 Select Report for Presentation screen" resides on the web server and the information needed to create the video is sent back to the Nurse Station.
- Select Reports Windows:: For a given CareGiver/ Care Recipient a list of the CR custom report templates (made previously). Default to highlighting the first report in the list.
- [View Prior Presentation] clears the current report selections and shows the reports and notes loaded for the prior presentation.
- View Reports Window: When a report is highlighted display the report/graph in the viewing window.
- Clinical Notes Window: Display the associated clinical notes for the selected report. These notes are only viewed by the nurse and are specific for the presentation/report only.
- Presentation Summary Window: Every presentation has a "Presentation Summary" report which the nurse prepares in outline format with bullet points. This Summary is saved as a report to be viewed in the presentation. This Summary Window allows the creation and editing of the summary.
- [RESET] resets the screen.
- [EXIT] exits the screen and returns to the Main NS Screen. Warn the user that exiting will loose all presentation prep work.
- [Next] saves the presentation info and goes to the next phase of screens: "2.1 Make Video" screen.
- [Back] returns to the prior screen.

Fig. 16

Screen Name: NS - 2.1 Make Video (resides on Nurse Station)



CareGiver: Eric Brown Care Recipient: Mary Smith		Presentation Date: mm/dd/yyyy Nurse: Ms. Kathy Kindess	
Clinical Notes for Selected Report (not included in presentation): Patient periodically is outside the recommended range.		02:31:00 Count Down 5 4 3 2 1 	
		Presentation Reports. Select the report for the presentation: Order [1] Presentation Summary [3] Report 1 [2] Report Weight/Temp [] Report Blood Pressure	
[REC] [PAUSE] [STOP]		[RESET] [PREVIEW]	
<div style="display: flex; justify-content: space-around; align-items: center;"> [<< BACK] NEXT >>] [EXIT] </div>			

Fig. 17

Overview:

1. Purpose of the screen is to make the video with the viewing of the reports
2. The video is recorded as Mirror Image.
3. [RECORD] does a count down from 5, before it starts recording
4. [PREVIEW] previews the last test recorded video. When first entering the screen, the [PREVIEW] is grayed-out and not accessible until a video has recorded.
5. The Presentation Summary Report and the previously selected presentation reports (from the prior screen) are displayed in the bottom right window as Thumb Nails. The ordering numbering of the reports remain blank. The ordering numbering indicates the order the reports were selected for presentation viewing. The system updates the ordering numbering after each report selection. When the Thumb Nail Report is selected the corresponding clinical notes and report are viewed by the nurse in the appropriate window. Note the Presentation Summary does not have any clinical notes associated with it.
6. The timing of the Thumb nail selection indicates when each report is to be displayed during the video viewing. Example:
 Step 1: Select Presentation Summary and report displayed in bottom left window. System places "1" as the order number.

Screen Name: NS - 3.1 Review Video and Flip Charts (resides on Nurse Station)



Presentation xxxx for <patient.fullname>	
	<div>00:31:00</div> 
<div>[Replay]</div>	
<div>[Deliver to Eric Brown] [REDO] [EXIT]</div>	

Fig. 18

Overview:

1. Purpose of the screen is view the presentation exactly as the CareGiver will view the presentation and submit the final presentation to be sent to the CareGiver.
2. [REPLAY] plays the presentation exactly how the presentation is to be viewed.
3. [Deliver to Eric Brown] delivers the presentation to the selected CareGiver through the web-server. Also, an alert is generated indicating a presentation has been sent and needs to be view. Lastly, an email is sent to CareGiver that a video has been prepared and needs to be review.
4. REDO returns the Nurse back to the NS-2.1 Make Video screen to redo the video.
5. EXIT exits the screen. Warn the user that exiting will loose all presentation prep work.