

BEST AVAILABLE COPY



10/29/2000

# Screen Name: 1.0 Home Page / Login

Screen Design Layout: Note: The screen design layout does not represent the Graphic Design nor the font selection and size.

TOP MENU	Home	CareHelper Guide	About Us	Contact Us	Site Map
<p>LEFT COLUMN MENU</p> <p> <input type="text"/> </p> <p>           If you are already a CareHelper User, welcome back!            Please enter your Caregiver User Name and Password.         </p> <p>           * Caregiver User Name <input type="text"/>             * Password <input type="password"/>   <input type="checkbox"/> (Save my password)         </p> <p> <input type="button"/> <input type="button"/> <input type="button"/> </p>					

\* Indicates required field entry

Legend: \* Indicates required field entry; <table,field name> defines field entry;  represents screen buttons;

Fig 1

FOOTER: 10/29/96

## Screen Name: 2.0 New Caregiver Registration

Screen Design Layout: Note: The screen design layout does not represent the Graphic Design nor the font selection and size.

TOP MENU Home Benefits of Membership Join Now TBD Contact Us

&lt; Error Messages Screen Area: Text to be displayed in red &gt;

## New Caregiver Registration

Just a few questions, and you are registered!

- \* Caregiver First Name
  - \* Caregiver Last Name
  - \* Caregiver Birth Date  [mm/dd/yyyy]
  - \* Caregiver Email Address
- Please check your Email Address! The confirmation of your registration will be sent to this address.

Choose your Caregiver User Name and Password. This will identify you in CareHelper.com

- \* Caregiver User Name
- ☐ (Save my User Name, so I don't have to remember it each time I login)
- Enter your password and re-type it to confirm. \* Indicates required field entry
- \* Password
- \* Confirm Password
- ☐ (Save my password, so I don't have to remember it each time I login)

Just one more thing.

As a registered Caregiver of CareHelper.com, you agree to abide by our Terms of Use. Please read the following Terms of Use before submitting your information.

\* Indicates required field entry

\* Indicates required field entry;

[table field name] defines field entry; [ ] represents screen buttons;



Fig 2



FOREF "FOZGTB60

## Screen Name: 3.0 Caregiver Page

Screen Design Layout: Note: The screen design layout does not represent the Graphic Design nor the font selection and size.

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us
<div> <div> LEFT COLUMN MENU (needs stable appcure)  ▼  &lt;caregiver.firstname&gt; Info  &lt;patient.nickname&gt; Info  Reports  Quick Glance  Personal Information  Contacts  Medical Information  Monitoring Devices  Medications  Supplies  Pharmacy / Store  ▼ </div> <div> <div> Welcome, &lt;caregiver.firstname&gt;,  to your personalized Home Page.  [Redacted] </div> <div> Select or Add a Care Recipient,  [Redacted] </div> <div> <div> Dad Joseph P. Smith, Sr. Inactive  &lt;patient.nickname&gt; &lt;patient.firstname&gt; &lt;patient.status&gt;  Add New Care Recipient [Redacted] </div> <div> <div> ALERTS </div> <div> &lt;patient.nickname&gt;  &lt;event.start&gt; &lt;event.description&gt;  &lt;event.start&gt; &lt;event.description&gt;  &lt;patient.nickname&gt;  &lt;event.start&gt; &lt;event.description&gt;  &lt;event.start&gt; &lt;event.description&gt; </div> </div> <div> <div> Shopping Manager </div> <div> &lt;caregiver.firstname&gt;'s List  &lt;shoplist.item_desc&gt;  &lt;shoplist.item_desc&gt;  &lt;patient.nickname&gt;'s List  &lt;shoplist.item_desc&gt;  &lt;patient.nickname&gt;'s List  &lt;shoplist.item_desc&gt;  &lt;shoplist.item_desc&gt; </div> </div> </div> </div> </div>					

Legend: \* Indicates required field entry.

&lt;table field name&gt; defines field entry.

[Redacted] represents screen buttons.

Fig 3

TOP LEFT TOP RIGHT

## Screen Name: 3.1 Caregiver Personal Information

Screen Layout Design:

TOP MENU Home

Janice's Organizer

Shopping Manager

CareHelper Guide

Contact Us

&lt; error messages screen area &gt;

## &lt;caregiver.nickname&gt;'s Personal Information

Please complete the following information and bring the information current:

*First Name	<caregiver.firstname>	*Last Name	<caregiver.lastname>
Nickname	<caregiver.nickname>	Birth Date	<caregiver.birthdate> (mm/dd/yyyy)
Gender	<caregiver.gender>		
Address	<caregiver.address1> <caregiver.address2>		
City	<caregiver.city>	Zip Code	<caregiver.zip>
State	<caregiver.state>		
*Day Phone	<caregiver.homephone>	Mobile Phone	<caregiver.mobile>
*Evening Phone	<caregiver.workphone>	Pager No.	<caregiver.pagerno>
Fax Phone	<caregiver.fax>	Email	<caregiver.email>
Primary Language	<caregiver.language1>	Secondary Language	<caregiver.language2>

Notes:

\* Indicates required field entry

[ &lt;caregiver.firstname&gt; ]

[ &lt;caregiver.lastname&gt; ]

[ &lt;caregiver.birthdate&gt; ]



Fig 4



102121 10291860

Screen Layout Design:

TOP MENU Home

Screen Name: 4.0.1 Add New Care Recipient

Janice's Organizer

Shopping Manager

CareHelper Guide

Contact Us

< error messages screen area >

## Add New Care Recipient

Welcome, <caregiver.firstname>.

### PROCESS STEPS

Personal Information

Emergency Contact Information

Medical Information

Monitor Devices

Medications

Supplies

Pharmacy / Store

Reports

Please complete the following information about your Care Recipient:

\* First Name

<patient.firstname>

\* Last Name

<patient.lastname>

\* Nickname

<patient.nickname> (eg. Mom, Dad, Nana, Patti, etc.)

\* Birth Date

<patient.birthdate> (mm/dd/yyyy)

\* Gender

<patient.gender>

\* Indicates required field entry

Fig. 5

Fig 5



Screen Layout Design:  
TOP MENU Home

Jarloe's Organizer

Shopping Manager

CareHelper Guide

Contact Us

< error messages screen area >

Screen Name: 4.1 Care Recipient Personal Information

<patient.nickname>'s Personal Information

PROCESS STEPS

Add New Care Recipient
Emergency Contact Information
Medical Information
Monitor Devices
Medications
Supplies
Pharmacy / Store
Reports

Please complete the following information about <patient.nickname> and bring the information current:

*First Name	<patient.firstname>	*Last Name	<patient.lastname>
*Nickname	<patient.nickname>	*Birth Date	<patient.birthdate> (mm/dd/yyyy)
*Gender	<patient.gender>	SSN	<patient.ssn> (123-45-6789)
Address	<patient.address1> <patient.address2>	Zip Code	<patient.zip>
City	<patient.city>	Evening Phone	<patient.phone2>
State	<patient.state>	Notification Email	<patient.email>
Day Phone	<patient.phone1>	Secondary Language	<patient.language2>
Mobile Phone	<patient.mobile>		
Pager No.	<patient.pagerno>		
Primary Language	<patient.language1>		

Notes:

\*Indicates required field entry

[ <patient.firstname> ]

[ <patient.lastname> ]

[ <patient.ssn> ]

Fig 6

T029T860

## Screen Name: 4.2 Care Recipient Emergency Contact Information

## Screen Layout Design:

TOP MENU

Home

Janice's Organizer

Shopping Manager

CareHelper Guide

Contact Us

&lt; error messages screen area &gt;

## &lt;patient.nickname&gt;'s Emergency Contact Information

## PROCESS STEPS

Please select an Emergency Contact and complete the following information:

## Select or Add an Emergency Contact,

Add New Care Recipient
Personal Information
Medical Information
Monitor Devices
Medications
Supplies
Pharmacy / Store
Reports

Primary Contact Order	Name	Relationship To Care Recipient
2)	Judith Conings	Neighbor
New Contact		

## \* Contact Order

&lt;contacts.order&gt;, in case of an emergency

## First Name

&lt;contacts.firstname&gt;

## Last Name

&lt;contacts.lastname&gt;

## Nickname

&lt;contacts.nickname&gt;

## Birth Date

&lt;contacts.birthdate&gt; (mm/dd/yyyy)

## \* Relationship

&lt;contacts.relationship&gt;

(i.e. Relationship to Care Recipient)

## \* Evening Phone

&lt;contacts.home&gt;

## \* Day Phone

&lt;contacts.work&gt;

## Mobile Phone

&lt;contacts.mobile&gt;

## Pager No.

&lt;contacts.pagerno&gt;

## Email

&lt;contacts.email&gt;

## Primary Language

&lt;contacts.language1&gt;

## Secondary Language

&lt;contacts.language2&gt;

Notes: (eg. Primary Contact has security code for home)

\* Indicates required field only





Fig 7



FOOTNOTES

### Screen Layout Design:

Screen Name: 4.3.3 Care Recipient Physician Information

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us																														
<p>&lt; error messages screen area &gt;</p>																																			
<p><b>PROCESS STEPS</b></p> <p>Add New Care Recipient</p> <p>Emergency Contact Information</p> <p>Medical History</p> <p>Monitor Devices</p> <p>Medications</p> <p>Supplies</p> <p>Pharmacy / Store</p> <p>Reports</p>																																			
<p><b>&lt;patient.nickname&gt;'s Medical - Physician Information</b></p> <p>Please select a Physician and complete the following information:</p>																																			
<table border="1"> <thead> <tr> <th>Name</th> <th>Select or Add a Physician, Specialty</th> <th>Status</th> </tr> </thead> <tbody> <tr> <td>Sally Jones, M.D.</td> <td>Geriatrics</td> <td>Inactive</td> </tr> <tr> <td colspan="3">Add New Physician</td> </tr> </tbody> </table>						Name	Select or Add a Physician, Specialty	Status	Sally Jones, M.D.	Geriatrics	Inactive	Add New Physician																							
Name	Select or Add a Physician, Specialty	Status																																	
Sally Jones, M.D.	Geriatrics	Inactive																																	
Add New Physician																																			
<table border="1"> <tbody> <tr> <td>Status</td> <td>&lt;physician.status&gt; &gt;</td> <td></td> </tr> <tr> <td>First Name</td> <td>&lt;physician.firstname&gt;</td> <td>Last Name &lt;physician.lastname&gt;</td> </tr> <tr> <td>Office Contact</td> <td>&lt;physician.contact&gt; &gt;</td> <td></td> </tr> <tr> <td>Office Phone</td> <td>&lt;physician.office&gt; &gt;</td> <td>Fax Number &lt;physician.fax&gt; &gt;</td> </tr> <tr> <td>Address</td> <td>&lt;physician.address1&gt;</td> <td></td> </tr> <tr> <td></td> <td>&lt;physician.address2&gt;</td> <td></td> </tr> <tr> <td>City</td> <td>&lt;physician.city&gt; &gt;</td> <td></td> </tr> <tr> <td>State</td> <td>&lt;physician.state&gt; &gt;</td> <td>Zip Code &lt;physician.zip&gt; &gt;</td> </tr> <tr> <td>Specialty 1</td> <td>&lt;physician.specialty1&gt; &gt;</td> <td>Other: &lt;physician.specialty1&gt; &gt;</td> </tr> <tr> <td>Notes:</td> <td colspan="2"></td> </tr> </tbody> </table>						Status	<physician.status> >		First Name	<physician.firstname>	Last Name <physician.lastname>	Office Contact	<physician.contact> >		Office Phone	<physician.office> >	Fax Number <physician.fax> >	Address	<physician.address1>			<physician.address2>		City	<physician.city> >		State	<physician.state> >	Zip Code <physician.zip> >	Specialty 1	<physician.specialty1> >	Other: <physician.specialty1> >	Notes:		
Status	<physician.status> >																																		
First Name	<physician.firstname>	Last Name <physician.lastname>																																	
Office Contact	<physician.contact> >																																		
Office Phone	<physician.office> >	Fax Number <physician.fax> >																																	
Address	<physician.address1>																																		
	<physician.address2>																																		
City	<physician.city> >																																		
State	<physician.state> >	Zip Code <physician.zip> >																																	
Specialty 1	<physician.specialty1> >	Other: <physician.specialty1> >																																	
Notes:																																			

Fig 8





## Screen Design Layout: Shopping Mgr

Screen Name: 3.2 Shopping Manager

TOP ENT Home

Janice's Organizer

Help/About/Screen

CareHelper Guide

Contact Us

&lt; error messages screen area &gt;

## Shopping Manager

## Suggestions and Reminders

<input type="button" value="SORT"/>	<input type="button" value="ADD ALL"/>	Item	Person	<input type="button" value="Delete"/>
11/23/1999	<item_description>	>	<patient,nickname>	<input type="checkbox"/>
11/19/1999	<item_description>	>	<patient,nickname>	<input type="checkbox"/>
11/18/1999	<item_description>	>	<patient,nickname>	<input type="checkbox"/>
11/18/1999	<item_description>	>	<patient,nickname>	<input type="checkbox"/>

## Prior Shopping Items

&lt;caregiver.firstname&gt;'s List

## Medications

&lt;item\_description&gt;

&lt;cg\_location&gt;

## Supplies

&lt;item\_description&gt;

&lt;cg\_location&gt;

## Other

&lt;item\_description&gt;

&lt;cg\_location&gt;

&lt;item\_description&gt;

&lt;cg\_location&gt;

&lt;patient,nickname&gt;'s List

## Medications

&lt;item\_description&gt;

&lt;cg\_location&gt;

&lt;item\_description&gt;

&lt;cg\_location&gt;

## Supplies

&lt;item\_description&gt;

&lt;cg\_location&gt;

\*\*\* CURRENT LIST \*\*\*

&lt;caregiver.nickname&gt;

Qty Item

001 <shoplist,item desc
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐

&lt;patient,nickname&gt;

Qty Item

001 <shoplist,item desc
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐
001 <shoplist,item desc	>	☐

&lt;patient,nickname&gt;

Qty Item

001 <shoplist,item desc
| 001 <shoplist,item desc | > | ☐ |
| 001 <shoplist,item desc | > | ☐ |

Fig 9



FOOTNOTES

Screen design Layout:

Screen Name: 4.7 Monitor Device Set Up

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us
< error messages screen area >					
<patient.nickname>'s Monitor Device Set Up					
Select or Add a Device,					
Device		Collection Method			
Pill Bottle		Auto	↓		
Blood Pressure		na			
Body Weigh		Manual			
Dose Pager		na			
New Device			↓		

Fig 10

### Screen design Layout:

**Screen Name: 4.7.5 Pill Bottle Device Set Up**

[illegible]


*Italic fields are display only (from medications set up screen)*

Fig 11

T02T2T T049T860

TOP MENU	Home	Janice's Organizer	Shopping Manager	CareHelper Guide	Contact Us
< error messages screen area >					
<b>&lt;patient.nickname&gt;'s Monitor Device Set Up</b>					
Select or Add Device,					
TBD		Collection Method			
		Device			
		Blood Pressure na			
		Blood Glucose Auto			
		New Device			
<u>Weight/Temp</u>					
Please select the date for <patient.nickname>'s Weight/Temp Set Up and complete the information:					
Weight/Temp <input type="text" value="&lt;Date&gt;"/> <input type="text" value="&lt;weight&gt;"/> [Add New Weight/Temp Measure]					
Please complete the following information:					
Measure Method Manual					
Date <input type="text" value="&lt;Date&gt;"/>					
Weight Measure <input type="text" value="&lt;wmeasure&gt;"/> Range: High <input type="text" value="&lt;weighthigh&gt;"/> Low <input type="text" value="&lt;weightlow&gt;"/>					
Temp Measure <input type="text" value="&lt;tmmeasure&gt;"/> Range: High <input type="text" value="&lt;temphigh&gt;"/> Low <input type="text" value="&lt;templow&gt;"/>					
		[Save and Continue]		[Finish Later] [Do Not Save]	

Fig 12



10216701-121201

Screen Name: NS - 1.1 Test Camera (resides on Nurse Station)

**Step by Step Instructions**

1. Adjust camera hardware
2. Preview video Adjust Volume, Contrast and Brightness
3. Preview Video
4. Save Setting




Fig. 14

**Overview:**

1. The purpose of this screen is to allow the Nurse to be able to view themselves and adjust the hardware and controls (Volume, Contrast and Brightness, only)
2. Test Camera Screen resides on the Nurse Station along with the Camera software.
3. The video is recorded as Mirror Image.
4. The camera is always on, see viewing the image starts immediately when first entering the screen.
5. [RECORD] does a count down from 5, before it starts recording
6. [PREVIEW] previews the last test recorded video.
7. Control settings show the current and recommended range of volume, contrast and brightness.
8. Above the Volume controls is a microphone recording feedback bar, which indicates the microphone recording level.
9. Note: The technical video settings that impact processing, quality, performance, compression and sending the video is to be accessed by a separate System Administrator screen only.

Screen Name: NS-1.2 Select CareGiver and Care Recipient (resides on web server)

Select CareGiver	Select Care Recipients
<caregiver.last>, <caregiver.first>	<Patient.initials>, <patient.first>, <patient.last>
Carroll, Pat	Dad
LeMore, Nancy	John
McCammon, Beth	Smith
Olters, Joan	
Thompson, Jane	
Victor, Andre	
Vace, Kenneth	
Young, Lee	
<div> <div>Continue</div> <div>EXIT</div> </div>	

Fig. 15

Overview:

1. Purpose of screen is to select who we are making the presentation for: CareGiver / Care Recipient
2. All of the CareGivers previously assigned to a nurse will be displayed in the CareGiver List. Default to highlighting the first CareGiver in the list.
3. When a CareGiver is highlighted, all of the Care Recipients associated with the CareGiver will be displayed in the Care Recipient window.
  - default highlighting the first Care Recipient listed.
4. This screen is based on the web-server, and accessed by the nurse through the Nurse Base Station. This interface is has seamless as possible to the nurse.
5. After selecting the CG, CR, and the Continue button is pressed, the selected information required to complete the Presentation is returned to the Nurse Station, to complete the video.
6. EXIT exits the screen and returns to the NS Main Screen.

Screen Name: NS - 1.3 Select Reports for Presentation (resides on web server)


CareGiver: Eric Brown Care Recipient: Mary Smith		Presentation Date: mm/dd/yyyy Nurse: Ms. Kathy Kindess	
Select Reports to include in Presentation: <input type="button" value="View Prior Presentation"/>		Presentation Summary Bullet Points: • Blood pressure in normal range • etc. •	
Include Reports <input type="checkbox"/> Report Template1 <input checked="" type="checkbox"/> Report Template2 <input checked="" type="checkbox"/> Report Template3 <input checked="" type="checkbox"/> Alerts	Description Report Template 1 Descr Report Template2 Report Template3	<input checked="" type="checkbox"/>	
View Selected Report:		Enter Clinical Notes for Selected Report: Patient periodically is outside the recommended range, and needs to be monitored closely during the holidays.	
			
<input type="button" value="BACK"/> <input type="button" value="NEXT"/> <input type="button" value="EXIT"/> <input type="button" value="RESET"/>			

Fig. 16




Overview:

1. Purpose of screen is to select the reports to be included in the presentation. In addition, clinical notes, that remain confidential to the clinician may be entered for each report
2. "1.3 Select Report for Presentation screen" resides on the web server and the information needed to create the video is sent back to the Nurse Station.
3. Select Reports Windows: For a given CareGiver/ Care Recipient a list of the CR custom report templates (made previously). Default to highlighting the first report in the list.
4. [View Prior Presentation] clears the current report selections and shows the reports and notes loaded for the prior presentation.
5. View Reports Window: When a report is highlighted display the report/graph in the viewing window.
6. Clinical Notes Window: Display the associated clinical notes for the selected report. These notes are only viewed by the nurse and are specific for the presentation/report only.
7. Presentation Summary Window: Every presentation has a "Presentation Summary" report which the nurse prepares in outline format with bullet points. This Summary is saved as a report to be viewed in the presentation. This Summary Window allows the creation and editing of the summary.
8. [RESET] resets the screen.
9. [EXIT] exits the screen and returns to the Main NS Screen. Warn the user that exiting will loose all presentation prep work.
10. [Next] saves the presentation info and goes to the next phase of screens: "2.1 Make Video" screen.
11. [Back] returns to the prior screen.



104121 10491860

Screen Name: NS - 2.1 Make Video (resides on Nurse Station)

Care Giver: Eric Brown Care Recipient: Mary Smith		Presentation Date: mm/dd/yyyy Nurse: Ms. Kathy Kindness	
Clinical Notes for Selected Report (not included in presentation): Patient periodically is outside the recommended range.		02:31:00	Count Down 5 4 3 2 1
			
		Presentation Reports. Select the report for the presentation: Order [1] Presentation Summary [3]  [2] Report Weight/Temp [ ] Report Blood Pressure	
[REC] [PAUSE] [STOP]		[RESET] [PREVIEW]	
[<< BACK] [NEXT >>] [EXIT]			

Overview:

1. Purpose of the screen is to make the video with the viewing of the reports
2. The video is recorded as Mirror Image.
3. [RECORD] does a count down from 5, before it starts recording
4. [PREVIEW] previews the last test recorded video. When first entering the screen, the [PREVIEW] is grayed-out and not accessible until a video has recorded.
5. The Presentation Summary Report and the previously selected presentation reports (from the prior screen) are displayed in the bottom right window as Thumb Nails. The ordering numbering of the reports remain blank. The ordering numbering indicates the order the reports were selected for presentation viewing. The system updates the ordering numbering after each report selection. When the Thumb Nail Report is selected the corresponding clinical notes and report are viewed by the nurse in the appropriate window. Note the Presentation Summary does not have any clinical notes associated with it.
6. The timing of the Thumb nail selection indicates when each report is to be displayed during the video viewing. Example:  
Step 1: Select Presentation Summary and report displayed in bottom left window. System places "1" as the order number.

Fig. 17

102121 10491860

Screen Name: NS - 3.1 Review Video and Flip Charts (resides on Nurse Station)



Presentation xxxx for <patient.fullname>	
	00:31:00
	
[Deliver to Eric Brown]	[REDO] [EXIT]

Fig. 18

Overview:

1. Purpose of the screen is view the presentation exactly as the CareGiver will view the presentation and submit the final presentation to be sent to the CareGiver.
2. [REPLAY] plays the presentation exactly how the presentation is to be viewed.
3. [Deliver to Eric Brown] delivers the presentation to the selected CareGiver through the web-server. Also, an alert is generated indicating a presentation has been sent and needs to be view. Lastly, an email is sent to CareGiver that a video has been prepared and needs to be review.
4. REDO returns the Nurse back to the NS-2.1 Make Video screen to redo the video.
5. EXIT exits the screen. Warn the user that exiting will loose all presentation prep work.

**This Page is Inserted by IFW Indexing and Scanning  
Operations and is not part of the Official Record**

## **BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☒ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINE(S) OR MARK(S) ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER: \_\_\_\_\_**

**IMAGES ARE BEST AVAILABLE COPY.**

**As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.**